

NUCLEAN - SURFACE DISINFECTION SPRAY
Ready-to-use spray for surface disinfection in medical environments



Product description:

NUCLEAN is a **ready-to-use surface disinfection spray** designed for medical, clinical and hospital environments. It is suitable for the **cleaning and disinfection of non-immersible medical device surfaces**, without rinsing (except for surfaces in contact with food).

The NUCLEAN spray combines **cleaning and disinfecting action** with fast drying and leaves no residue. Its **hypoallergenic formulation** makes it suitable for frequent use in healthcare settings.

Main features

- Ready-to-use spray
- For non-immersible medical device surfaces
- No rinsing required (except food-contact surfaces)
- Fast drying, streak-free
- Hypoallergenic formulation
- Complete disinfectant activity in 5 minutes

Packaging

- **500 ml spray bottle**
- Can be used with non-woven wiping cloths (e.g. bionet wipes 40 x 38 cm, ref. 712)

Microbiological properties - NUCLEAN Spray

Bacteria

- EN 13727: active in 1 minute
- EN 1276: active in 5 minutes
- EN 13697: active in 5 minutes
- EN 14561: active in 5 minutes

Yeasts / Moulds

- EN 13624: active in 1 minute
- EN 1650: active in 5 minutes
- EN 13697: active in 5 minutes
- EN 14562: active in 5 minutes

Viruses

- EN 14476: Norovirus, HIV-1 – active in 30 seconds
- EN 14476: Hepatitis B, Herpes virus, Influenza A/H1N1 – active in 1 minute
- EN 14476: Hepatitis C, Rotavirus, Adenovirus – active in 5 minutes

Spores

- EN 13704: *Clostridium difficile* – active in 15 minutes
- NF T 72-281: *Bacillus subtilis* – active in 5 minutes

Directions for use - NUCLEAN Spray

- 1Pre-impregnate a clean non-woven wipe with the solution.
- 2Thoroughly clean the surface.
- 3Leave to act for **5 minutes** (15 minutes for *Clostridium difficile*).
- 4Do not rinse, except for food-contact surfaces.
- 5Use a new wipe for each piece of equipment.

Reference: EC500