

STERANIOS 2% - COLD DISINFECTANT FOR MEDICAL DEVICES High-level cold disinfection solution for medical and endoscopic devices



Product description:

STERANIOS 2% is a **ready-to-use cold disinfection solution based on glutaraldehyde**, designed for the **high-level disinfection of heat-sensitive medical devices**. It is used by immersion in a dedicated soaking tray for surgical, medical and endoscopic equipment.

The solution must be used **in combination with a pre-cleaning disinfectant detergent**. The soaking bath can be **kept for up to 7 days**, provided that the glutaraldehyde concentration is regularly checked using reactive test strips.

- **High-level disinfection by immersion**
- **Ready-to-use solution**
- **Suitable for heat-sensitive medical devices**
- **Use in a dedicated soaking tray**

Microbiological activity

- **Bactericidal** (EN 13727, EN 14561, Helicobacter pylori): active in 10 minutes
- **Mycobactericidal** (EN 14348, EN 14563, Mycobacterium massiliense): active in 10 minutes
- **Fungicidal** (EN 13624, EN 14562): active in 10 minutes
- **Virucidal** (EN 14476): active in 10 minutes
- **Sporicidal**: demonstrated after 1 hour at 20 °C (NF T 72-230/231, EN 13704)

Directions for use

- 1Pre-disinfection:** clean the medical device using a suitable pre-disinfectant detergent, then rinse thoroughly. For endoscopic equipment: brush and irrigate.
- 2Disinfection:** pour STERANIOS 2% into a dedicated soaking tray and cover. Contact time: **10 minutes**.
- 3Rinsing:** rinse thoroughly with sterile or filtered water (0.2 µm), inside and outside for endoscopes.
- 4Drying and storage:** wipe with a single-use sterile cloth and store aseptically until next use.

Disinfection level: Level 3 (high)

Reference: AN382034