

SEPTIWIPES - DISINFECTION WIPES LEVEL 1
Alcohol-free wipes for low-level disinfection of medical devices



Product description:

SEPTIWIPES are **alcohol-free impregnated wipes** designed for the **cleaning and low-level disinfection** of **non-invasive reusable medical devices**. They are particularly suitable for **alcohol-sensitive equipment**, especially **superficial ultrasound probes**.

Ready to use, SEPTIWIPES integrate easily into wipe-based disinfection protocols in medical practices, imaging departments, clinics and healthcare facilities. Their formulation preserves sensitive materials while ensuring microbiological efficacy in compliance with applicable standards.

Key features

- Alcohol-free disinfectant wipes
- For low-level disinfection of non-invasive medical devices
- Recommended for superficial ultrasound probes
- Wipe size: **140 × 200 mm**
- **Neutral pH (pH 7)**
- Ready to use, **no rinsing required**

Microbiological activity

Bacteria & yeasts

- EN 16615, EN 1040, EN 13727, EN 1275, EN 13624 (dirty conditions): **effective in 1 minute**

Viruses

- EN 14476 (clean conditions) & RKI/DVV protocol: **effective in 30 seconds**
- EN 14476 (dirty conditions): **effective in 2 minutes**

Ultrasound probe compatibility

SEPTIWIPES are certified for ultrasound probe disinfection by the following manufacturers: **GE, Siemens, Philips, Canon Medical, Shimadzu, SonoAce, Aloka, Esaote, Hitachi.**

Regulatory information

- Medical device **Class IIa** (EU Regulation 2017/745)
- **Non-sterile**, professional use only
- Alcohol-free, latex-free, phthalate-free
- Shelf life after opening: **3 months** (lid tightly closed)
- Storage conditions: **+5 °C to +25 °C**

Packaging

- Dispenser box of **100 wipes**
- Reference: **SP2000**

Directions for use

- 1 Remove one wipe from the dispenser.
- 2 Thoroughly wipe the surface or medical device.
- 3 Allow the product to act for the required contact time.
- 4 Do not rinse.
- 5 Use a new wipe for each device.

Disinfection level: Level 1 (low)

Reference: SP2000