

RING PESSARY - UTERINE PROLAPSE STAGE 1 AND 2 Regular or thick ring pessary.



Product description:

The ring pessary is made of high quality silicone. It can be stored at temperatures of 1 to 30 ° C away from UV rays without contact with reactive media, gas, ozone or mineral oil. All models have been adjusted by a compression test so that the same craftsmanship is required to bend the device regardless of size.

Application

The annular pessary is elastic and can be easily changed by the patient herself. The circumference of the ring is relatively thin. Preferably, the mild forms of genital prolapse and incontinence are treated by the ring pessary. An intact pelvic floor (even if reduced) is required for this pessary. In addition, it can prevent the development of stress incontinence by its suburethral support effect.

Note

The maximum period of uninterrupted application is 30 days. It is recommended that the patient changes the pessary herself, that is, she should remove it every evening and insert it in the morning. The patient should be trained in the application / change by a healthcare professional, whereby additional treatments (eg creams) should be discussed. The product should not be used by other patients. The pessary can be cleaned under running water without using disinfectant.

Size: Regular, Thick

Pessary size: Size 1 - Ø50 mm, Size 2 - Ø55 mm, Size 3 - Ø60 mm, Size 4 - Ø65 mm, Size 5 - Ø70 mm, Size 6 - Ø75 mm, Size 7 - Ø80 mm, Size 8 - Ø85 mm, Size 9 - Ø90 mm, Size 10 - Ø95 mm, Size 11 - Ø100 mm



MEDICAL EQUIPMENT, ACCESSORIES AND CONSUMABLES FOR MEDICAL IMAGING.

Reference: DR0