

**DR. ARABIN RING PESSARY – UTERINE PROLAPSE STAGE 1 AND 2**  
**Silicone vaginal ring pessary for mild genital prolapse and incontinence**



**Product description:**

The Dr. Arabin ring pessary is a medical device made of **high-quality silicone**. It is intended for patients with **mild forms of genital prolapse** and certain cases of **urinary incontinence**, provided that the pelvic floor is strong enough to retain the device.

Thanks to its elasticity, the ring pessary can be **easily inserted and removed by the patient herself**, after appropriate training by a healthcare professional. All models are adjusted by **compression testing** to ensure similar handling regardless of the size.

The pessary is designed for **single-patient use** and must not be used by other patients. It can be **cleaned under running water** without the use of disinfectants.

**Indication / Intended use**

The ring pessary is indicated for patients with mild genital prolapse and incontinence. The pelvic floor must be sufficiently strong to maintain the pessary in position.

The maximum uninterrupted application period is **30 days** (medical device **Class IIa**). It is recommended that the patient removes the pessary in the evening and reinserts it in the morning, after being instructed by a healthcare professional.

**Storage conditions**

The pessary can be stored at temperatures between **0 and 50 °C**, protected from UV radiation and without contact with reactive media, gas, ozone or mineral oil.

**Available sizes**

The Dr. Arabin ring pessary is available in the following diameters:

50 mm, 55 mm, 60 mm, 65 mm, 70 mm, 75 mm, 80 mm, 85 mm, 90 mm, 95 mm, 100 mm.

Size: Regular

Pessary size: Size 1 - Ø50 mm, Size 2 - Ø55 mm, Size 3 - Ø60 mm, Size 4 - Ø65 mm, Size 5 - Ø70 mm, Size 6 - Ø75 mm, Size 7 - Ø80 mm, Size 8 - Ø85 mm, Size 9 - Ø90 mm, Size 10 - Ø95 mm, Size 11 - Ø100 mm

Reference: -