

HEAD PROTECTION - PRONE HEADREST Head gel positioner for operating table.



Product description:

Thanks to its particular ergonomic design, it can be used as a neck and head support during all types of surgery, with an easy access for intubation.

Materials

It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm³.

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between -18 and 50 degrees
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

Package

1 Piece.

Size of gel positioning wedges 1: Choose an option, 28 x 23 x 14,5 cm, 23 x 19 x 13 cm

Reference: -