

HEAD PROTECTION - ADULT/KID PRONE HEADREST Adult/kid Prone Headrest



Product description:

The gel positioners are devices that are used to reach an optimal patient position on the operating table. Many differentes models are available in various shapes and sizes; all fully compatible with human tissues. Some of the main features are: biocompatibility, comfort and extraordinary mechanical resistance.

Description

The PRONE HEADREST is made of an amorphous polyurethane gel with a density of 0.99 gr/cm3. Thanks to its particular ergonomic design, it can be used as a neck and head support during all types of surgery, with an easy access for intubation.

Materials

Amorphous Polyurethane Gel with a Density of 0.99gr/cm3

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room
- Avoid immersing the gel positioners in the detergent and / or disinfectantfor a long time

Storage and Disposal

- Keep the positioners, if possible, on a level surface
- Avoid placing heavy or sharp objects on the positioners
- Keep the gel temperature between -18 and 50 degrees
- Keep away from sources of light, moisture and heat
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

<u>Package</u>

1 piece.



Size of gel positioning wedges 1: $25 \times 20 \times 10$ cm - Adult, $20 \times 16 \times 11$ cm - Kid, Choose an option

Reference: CDP0003