


HEAD PROTECTION PAD - HORSESHOE WITH CENTRAL MEMBRANE Gel head and neck positioning wedges for operating table



Product description:

It can be used as a neck, head and face support during all types of surgery, with an easy access for intubation.

Materials

It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm³.

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between 50 ° and -18 °.
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

Package

1 Piece.

Size of gel positioning wedges 1: Choose an option, 9 x 4 x 2 cm, 14 x 6 x 3,5 cm, 20 x 7,5 x 4,5 cm, 20 x 7,5 x 7,5 cm

Reference: CDP0006