

SPOR'ACTIV - LIQUID DISINFECTANT LEVEL 3 Ultrasound probe disinfection liquid



Product description:

More information about disinfection levels [here](#).

High level disinfection in 3 minutes !!!

- Use :
 - Step 1: Pre-treatment with the use of a Clean'Up wipe
 - Step 2: High level disinfection by taking another Clean'Up wipe and applying one pressure of Spor'Activ to it.
- Clean'Up: Alcohol-free wipes
- Spor'Activ: Stabilized peracetic acid solution

Characteristics

- 1 liter container

Disinfection level: Level 3 (high)

Reference: -