


**SOAKING TANK FOR DISINFECTION OF PROBES**  
**Immersion disinfection tray for endocavitary probes**



**Product description:**

This **soaking tray** is designed for the **disinfection of endocavitary probes** by immersion in an appropriate **disinfectant solution**. It allows effective disinfection of probes used in medical imaging, particularly in gynecology and urology.

When used with the **SEPTICIDE** disinfectant solution, this tray enables **high-level disinfection in 10 minutes**, in accordance with immersion disinfection protocols.

**Features**

- Cylindrical tray with lid
- Immersion disinfection for endocavitary probes
- High-level disinfection in **10 minutes** with the appropriate solution
- Integrated hooks to secure cables and connectors
- Can be used **standing or wall-mounted**
- Usable volume: **2.5 liters**
- Replacement trays available

**Compatibility**

This tray is intended for use with the SEPTICIDE **disinfectant solution**, in accordance with the manufacturer's instructions and the facility's internal protocols.

**Use**

- Fill the tray with the appropriate disinfectant solution

- Immerse the endocavitary probe according to the protocol in force
- Observe the required contact time for disinfection
- Remove the probe and proceed with the next steps according to good practice

Config bac trempage: Disinfection tray with support, Set of 3 replacement trays

Reference: -