

**SEPTISPRAY - DISINFECTION SPRAY LEVEL 1**  
**For surfaces and medical devices cleaning and disinfection.**



**Product description:**

More information about disinfection levels here.

Apply by spraying on surfaces (benches, trolleys, furniture), equipment (monitors) and non-submersible, invasive and non-invasive medical devices (ultrasound probes, endocardial leads and sheathed endoscopes, Hemodialysis generators, dental unit).

Rapid action: Activity's spectrum achieved in just 5 minutes.

- Eliminate 99, 9999% of bacteria.
- Effective on resistant Clostridium spores.
- Material's compliance : alcohol free, neutral pH
- User's compliance : fragrance free, allergen free

**Formulation :**

Liquid preparation, TP2/TP4 operational.

DIDECYL DIMETHYLAMMONIUM CHLORIDE 7173-51-5 0.25 g/kg. ALKYL BENZYL

DIMETHYLAMMONIUM CHLORIDE 68424- 85-1 0.95 g/kg

**Features :**

Physical condition: Liquid

pH : 7,0 +/- 0,5

Perfume : None

Color : Colorless to yellow

**Precautions for use :**

Use the biocide product with care. Before use, read the label and the information relative to the product. MD Classe

IIb.

**Characteristics**

- Ultrasound probe cleaner

Disinfection level: Level 1 (low)

Reference: PROBE51