

SONY PRINTER UP-X898MD B/W A6 Compact medical grade black and white digital printer



Product description:

High quality printing:

Advanced thermal printing technology enables photo-quality black & white prints with a high resolution up to approximately 325 dpi and grey levels up to 256 steps.

A6 prints in less than two seconds:

While operating in high speed mode, these printers can produce A6-size prints in approximately 1.9 seconds (when used with UPP-110HG or UPP-110S print media).

Accepts both digital and analog video input signals:

A choice of analog composite video*3 and digital inputs allows flexible use with a wide range of modality devices.

USB flash drive storage*4:

The UP-X898MD can automatically store printed images from both USB and Video inputs onto a connected USB flash drive.

Compact, space-saving design:

Measuring 88mm (H) x 154 mm (W) x 240 mm (D), these printers have the same space-saving design as the predecessors for easy, flexible integration into medical carts.

Enhanced operability

Building on the success of the popular UP-897MD and UP-D897, these printers feature an enlarged LCD panel that displays printer settings, selected paper type and advisory messages. An intuitive joystick control allows speedy menu navigation, while the new Auto Lock mode prevents unintentional changes to print settings.

Hard copy reference of print settings:

For convenient reference, selected settings (Brightness, Contrast, Gamma and Sharpness) can optionally be printed in a blank space on the print paper.

Uses widely available print media:

Sony offers a choice of three high quality, widely available thermal print media types to meet the needs of different applications. Choose from the UPP-110HD, UPP-110S, or the UPP-110HG High Glossy Printing Paper that produces glossy, long-lasting near-photo quality prints.

Wide PC compatibility:

For maximum flexibility, a wide range of PC drivers are provided for Windows XP, Vista, 7 and 8 (32/64bit).

Quick guide print:

Quick instruction guides for different print operations can be selected via the front control panel. They can be printed out for easy reference without needing to refer to the user manual

Color: Black, White

Reference: -