

STERILET INSTALLATION SET Single Use Sterile



Product description:

To improve your daily practice, Ima-x offers you an IUD insertion set consisting of:

- The straight Pozzi forceps allow the cervix to be grasped and the cervix to be "traction-fixed" to allow the uterine flexion angle to be adjusted to 180°.
- The angled Chéron forceps clean the uterus and the work area.
- The lateral or distal thread cutting scissors are used to cut the threads of the intrauterine device.
- The hystrometer is used to assess the depth of the uterine cavity.

Characteristics

- 1 Hystrometer
- 1 Pozzi forceps
- 1 wire-cutting scissors
- 1 Cheron forceps

Reference: -



MEDICAL EQUIPMENT, ACCESSORIES AND CONSUMABLES FOR MEDICAL IMAGING.