

**OVAL PATIENT POSITIONING PAD
Gel Positioner for Operating Table.**



Product description:

This extremely versatile positioners are good for any procedure. Beneficial for tracheotomy, thyroidectomy, tonsil and adenoid cases, cataract or other ophthalmic procedures.

Materials

It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm³.

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between 50 ° and -18 °.
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

Package

1 Piece.

Size of gel positioning wedges 1: Choose an option, 16 x 8 x 4.5 cm, 33 x 7 x 5 cm, 33 x 12 x 6 cm, 36 x 12.5 x 7 cm, 25 x 12.5 x 5 cm, 38 x 12.5 x 6.3 cm, 50 x 12.5 x 7 cm

Reference: -