

GEL POSITIONING PAD - SHOULDER
Amorphous polyurethane gel pads



Product description:

The gel positioners are devices that are used to reach an optimal patient position on the operating table. Many different models are available in various shapes and sizes; all fully compatible with human tissues. Some of the main features are: biocompatibility, comfort and extraordinary mechanical resistance.

Materials

Amorphous PolyurethaneGel with a Density of 0.99 gr/cm³

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time

Storage and Disposal

- Keep the positioners, if possible, on a level surface
- Avoid placing heavy or sharp objects on the positioners
- Keep the gel temperature between 50 ° and -18 ° Keep away from sources of light, moisture and heat
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

Package

1 Piece

Size of gel positioning wedges 1: Choose an option, 30 x 22 x 5/1 cm, 50 x 12 x 10/4

Reference: -