

## RELTON-HALL FRAME POSITIONING PADS SET Relton-Hall Gel Positioning Wedges



### Product description:

This 4-piece set is designed to fit snugly over the four posts of the Relton-Hall Frame. Patients in the prone position on the Relton-Hall Frame, have extreme pressure localized at the four contact areas of the posts. The Relton-Hall Pad Set provides pressure protection, reducing the chance of pressure related tissue/nerve damage to the patient, regardless of surgery duration.

### Materials

It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm<sup>3</sup>.

### Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

### Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between 50 ° and -18 °.
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

### Package

4 Pieces.

Size of gel positioning wedges 1: 19 x 17.5 x 10.5

Reference: -