

GEL HEEL SUPPORT Gel positioning pad- 2 pieces



Product description:

It should be used for any procedure that is lengthy and involves pressure on the patient's heel area.

Materials

It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm³.

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between 50 ° and -18 °.
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

Size of gel positioning wedges 1: Choose an option, 20 x 11 x 7 - Adult, 9.7 x 5.3 x 3.5 - Pediatric, 7.8 x 4.3 x 2.8 - Neonatal

Reference: -



MEDICAL EQUIPMENT, ACCESSORIES AND CONSUMABLES FOR MEDICAL IMAGING.