

**POSITIONING MATTRESS - PERINEAL SHAPE**  
**Perineal mattress pad for operating table in amorphous polyurethane gel**



**Product description:**

Thanks to the various sizes, they are designed to cover the operating table during any surgical procedure. The perineal cut allows a complete discharge at the sacral level. They are suitable for all those patients who need to protect this body area.

Materials

It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm<sup>3</sup>.

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between 50 ° and -18 °.
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

Package

1 Piece.

Reference: -