

PROFILED POSITIONING WEDGE - ULNAR NERVE Gel wedge for arm positioning



Product description:

It is anatomical and it can be used as arm support during all types of surgery. It is specially designed for the protection of the ulnar nerve.

Dimensions: 52 x 15 x 5.5 cm

Materials

It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm³.

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between -18 and 50 degrees.
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

Package

1 Piece.

Reference: -