

Required Information

When Requesting an offer for a Protective Curtain Model VH194-NV

0	Area of Use:	Indoor X-ray protection for windows (Not suitable for outdoor use or as protection for doors or entrances.)	
1.0	Curtain Size: (If known)	cm x cm (total height including track x complete length/track length) (Fill out only after taking steps 2.1 to 4.0 into consideration . Steps 2.0, 5.0 to 6.2, and 7.0 still need to be completed!)	
2.0	MountedTo:	☐ Ceiling or ☐ Wall	
2.1	MountingVariation:	☐ Inside window recess ☐ Outside window recess	
2.2	FurtherInformation:	cm Recess height cm Recess width	
2.3	Spacing Left: Right:	cm Recess wall to window cm Room wall to window cm Window to recess wall cm Window to room wall	
		(With enough space, the curtain will not obstruct the window in a parked state.)	
3.0	Window Size: (Including frame)	cm x cm (H x W) cm x cm (H x W)	
4.0	Window Position:	cm Upper edge of window to upper edge of recess cm Lower edge of window to lower edge of recess cm Lower edge of window to lower edge of recess cm Lower edge of window to floor	
5.0	Curtain Pb-Value: (Maximum curtain length)	□0.50 mm □1.00 mm □1.50 mm □2.00 mm (6.80 m) (4.08 m) (2.82 m) (2.04 m)	IMA-X s.a. Equipment, accessorie and consumables for
6.1	Panel Park Position:	Right or Left (When possible, the park position should be outside of the window.)	medical imaging Adress:
6.2	Position of Movement Cords:	Right or Left	Am Pëtz 17 L-9579 Weidingen Luxemburg
7.0	Additional Information: (Optional)	Please give the position and measurements of any special structural conditions, radiator, protruding window bank, etc.	Phone: +352 95 75 52 +352 95 97 45
	Important Information	on	e-Mail: info@ima-x.com
	_	For safety, there should be an overlap of 10 cm on all sides	Web: www.ima-x.com
		2-3 cm (min. 1 cm) from the bottom of the panels to the floor nat the panels can properly move and turn.	