

# 3-LUMENS MALLEABLE HYSTERO CATHETER WITH BALLOON PM-CARE® 3V-BALLOON





## **Product description:**

#### **Description:**

The PM-CARE® 3V-BALLOON 3-lumens malleable hystero catheter with balloon is used to perform hysterography and hysterosonography. This malleable catheter can be bent and shaped to adapt to the patient's anatomy.

### Advantages: • Easy and simple manipulation.

- Sealing guaranteed by the silicone balloon.
- Standard connector for connecting a syringe for product injection.
- Atraumatic insertion.

Catheter diameter: 2.66 mm / 8 Fr

Inflated Balloon diameter : 13 mm

Balloon insufflation max gas : 2.5 cc

Length: 230 mm

#### **About this Medical Device**

1st release date: 2011

Single use: Yes





Class: Is

Sterile: Yes

Sterilisation method: Ethylene oxide sterilisation

Contains latex: No

Lifetime: 5 years

Certification: CE

Packaging: Box of 5 units in individual sterile wrapping

Storage: In its original box in a cool and dry place

Reference: PS3310733