

KRYSTAL – STERILE POLYURETHANE PROBE COVER WITH INTEGRATED ADHESIVE GEL FOR SUPERFICIAL PROBES
Sterile superficial ultrasound probe cover with integrated adhesive gel



Product description:

Krystal polyurethane probe covers with integrated adhesive gel are designed to provide **long-lasting and effective protection for superficial ultrasound probes** during ultrasound-guided procedures. They help limit exposure of the probe to **bodily contaminants** and facilitate **cleaning and disinfection** of the equipment after use.

The **integrated sterile adhesive gel** at the tip of the cover allows precise positioning on the probe without the need to apply gel inside the cover. This design improves probe contact while simplifying exam preparation.

The covers are supplied as **ready-to-use sterile kits**, including all accessories required for the procedure, ensuring **optimal comfort for medical staff** and compliance with aseptic protocols.

Key benefits

- **Time-saving and practical:** no need to apply gel inside the cover, reducing gel residue on the probe
- **Artifact-free imaging:** adhesive gel tip adapts to the probe shape, preventing folds and air bubbles
- **Latex-free material:** polyurethane does not cause irritation or intolerance

Each box contains **20 individually packaged sterile kits**, ready for use.

Available references

Reference	Description	Dimensions	Packaging
YAPC1291	Cover with adhesive gel, elastic bands and adhesive strips	+13 x 61 cm	20 kits / box
YAPC1290	Cover with adhesive gel, elastic bands and adhesive strips	+15 x 120 cm	20 kits / box

YAPC1292

Cover with adhesive gel +15 x 144 cm
gel, elastic bands and
adhesive strips

20 kits / box

Contents of the sterile kits

- 1 sterile Krystal probe cover with integrated adhesive gel for superficial probes
- 1 EDM sterile gel sachet, 20 g
- 2 fixation elastic bands
- 1 sterile drape

Choosing the appropriate length

Covers **shorter than 60 cm** are suitable for short outpatient procedures such as injections, infiltrations and vascular access. Covers measuring **60 to 125 cm** are recommended for biopsy and anesthesia procedures. Covers **longer than 125 cm** are intended for use in operating rooms.

The use of a probe cover **does not replace proper disinfection in accordance with current good practice guidelines.**

Size: 13 x 61 cm, 15 x 120 cm, 15 x 144 cm

Ultrasound probe covers - Material: Polyurethane

Gel Protections sondes: With gel

Reference: -