


## OVERHEAD CEILING MOUNTED SHIELD 350 Eyes and thyroid radiation shields model 350



### Product description:

- Material: Lead glass, with minimum 2.2 mm lead equivalent.
- Fitted flexible lower curtain: lead equivalent 0.5 mm.

40 cm high x 50 cm wide shield with a suspension arm and ceiling column

40 cm high x 50 cm wide shield with a suspension arm, column and 250 cm ceiling track

30 cm high x 40 cm wide shield with a suspension arm and column

30 cm high x 40 cm wide shield with suspension arm, column and 250 cm ceiling track

Also available with a led lamp :

Stationary system: with a fixed ceiling column

40 x 50 cm shield, static column and lamp LED 130F

40 x 50 cm shield, static column and lamp LED 150F

40 x 50 cm shield, static column and lamp LED 300 DF

Movable solution: with ceiling column and 250 cm HD aluminium ceiling track

40 x 50 cm shield, column, ceiling track and lamp LED 130F

40 x 50 cm shield, column, ceiling track and lamp LED 150F

40 x 50 cm shield, column, ceiling track and lamp LED 300 DF



MEDICAL EQUIPMENT, ACCESSORIES AND CONSUMABLES FOR MEDICAL IMAGING.

Reference: -