

LEADED APRON 527

Lead apron closed in front by one panel that fasten above the chest.



Product description:

Lead apron closed in front by one panel that fasten above the chest. Shoulder padding for optimal wearing comfort. An inner pocket is stitched on the liner that is specially designed for storing a dosimeter during interventions. This model can be supplemented by a wide support belt to better distribute the weight of the apron.

Options available:

- Wide belt
- Strap belt
- Two-tone apron
- Outside breast pocket
- Non-removable sleeve
- Embroidery (30 characters)
- Embroidery on velcro (15 characters)
- Embroidery logo
- RFID Tag

Size: XS, S, M, L, XL, XXL

PA Color: 9. Grey, 10. Neon yellow, 11. Neon orange, 12. Green, 13. Royal blue, 14. Midnight blue, 15. Plum

Front Pb equivalent: 0.25 mm Pb, 0.35 mm Pb, 0.50 mm Pb

Rear Pb equivalent: 0.25 mm Pb, 0.35 mm Pb

Closing - x-ray protection clothing: Velcro, Buckle

Material - x-ray protection clothing: PA - Polyamide, PVC

Protective range: S - Light, U - Ultra-light

PVC color: 1. Purple, 2. Lilac, 3. Indian rose, 6. France Blue, 7. Sky Blue, 8. Green, 5. Navy blue, 4. Red

Reference: VDR0005