

LEADED HALF-CHASUBLE - BLOC
Radiation protection half-chasuble



Product description:

This frontal lead apron is made up at the front of a panel which is held to the body by two panels which cross at the back.

The braces are not stitched under the arm, so that this model, specially designed for surgeons, can be unfastened and removed without touching the sterile gown.

Once the fastenings have been undone, the apron slips off by itself.

The shoulders are padded for improved comfort, an inside pocket sewn into the lining has been specially designed to hold your dosimeter.

With our system of loops, the apron can be hung up easily on our hangers and specially made supports.

Options:

- P-Thy incorporated Lead Equivalent 0,35mm
- P-Thy incorporated Lead Equivalent 0,50mm
- Velcro fastening
- Two-tone apron
- Outside breast pocket
- Embroidery (30 characters)
- Embroidery on velcro (15 characters)
- Embroidery logo
- RFID Tag

Size: XS, S, M, L, XL, XXL

Lead equivalent: 0.35 mm Pb, 0.50 mm Pb, 0.25 mm Pb

PA Color: 19. Grey, 13. Green, 16. Petrol blue, 17. Black, 18. Plum, Lavender, 15. Red

Material - x-ray protection clothing: PA - Polyamide, PVC

Protective range: S - Light, U - Ultra-light

PVC color: 8. Purple, 9. Lilac, 12. Red, 5. Navy blue, 7. France Blue, 4. Khaki, 3. Taupe, 1. Grey, 2. Charcoal, 6. Teal, 10. Raspberry, 11. Dusty pink

Reference: -